



Sri Dharmasthala Manjunatheshwara College
(Autonomous), Ujire-574 240, Dakshina Kannada, Karnataka State

2.3.2

ICT ENABLED TEACHING SUPPORTIVE DOCUMENT

LAB MANUAL

DEPT OF CHEMISTRY (UG)

PRACTICAL MANUAL

B.Sc - II Semester

2019-20

Name:.....

Seat No:.....



Department of Chemistry
SDM College
Ujire

INDEX

Sl.No	Name of the experiment	Page no
1	Calibration of pipette	2-3
2	Estimation of sodium hydroxide	4-9
3	Estimation of hydrochloric acid	10-15
4	Estimation of Mohr's salt	16-21
5	Estimation of ferric chloride	22-27
6	Estimation of ferrous and ferric ions in a mixture	28-33
7	Estimation of copper	34-39
8	Estimation of oxalic acid and sulphuric acid in a mixture	40-45
9	Estimation of calcium in chalk	46-49
10	Determination of hardness of water	50-57
11	Estimation of glucose	58-63
12	Estimation of vitamin C	64-69
13	Determination of strength of sodium hydroxide by two burette titration	70-71

OBSERVATIONS AND CALCULATIONS

Laboratory temperature	=	$^{\circ}\text{C}$
Density of water at laboratory temperature	=	gcm^{-3}
Mass of empty beaker	= w_1 =	g
Mass of beaker + water	= w_2 =	g
Mass of water	= $(w_2 - w_1) \text{ g}$	

=

=

g

$(w_2 - w_1)$

Volume of water pipette out to the beaker = $\frac{(w_2 - w_1)}{\text{Density}} \text{ cm}^3$

=

=

cm^3

Expt.No:

Date:

Calibration of pipette

Aim: To calibrate the pipette

Procedure:

The pipette provided is cleaned using chromic acid. This is done by keeping the pipette in chromic acid for about an hour. The pipette is then thoroughly washed with water and rinsed with distilled water.

Then the pipette is filled with distilled water which is at room temperature. Water is run out until the lower meniscus is exactly on the mark. The drop adhering to the tip of the pipette is removed by the surface of water contained in a beaker in contact with the tip.

The pipette is then allowed to discharge into a clean and dry beaker which is previously weighed. The tip of the pipette should be in contact with the side of the beaker. The pipette is allowed to drain for 15 seconds after the outflow has stopped. Now the pipette is removed from the beaker. Care is taken that no drop of water is adhering to the outside of the pipette.

The beaker containing the water collected from the pipette is weighed again. The room temperature is noted.

Result:

The calibrated volume of the pipette = cm^3

OBSERVATIONS AND CALCULATIONS

(a) Preparation of standard sodium carbonate solution

Mass of weighing bottle + sodium carbonate = w_1 = g

Mass of weighing bottle after transferring = w_2 = g

Mass of anhydrous sodium carbonate transferred = $(w_1 - w_2)$ g

=

=

g

Strength of sodium carbonate solution

$$= \frac{(w_1 - w_2) \times 4}{53}$$

=

=

N

(b) Standardization of hydrochloric acid solution

Titration of 25cm^3 of standard sodium carbonate solution against hydrochloric acid solution

Indicator : Methyl orange

End point : Yellow to orange

Trial no.	I	II	III.
Burette reading			
Final reading			
Initial reading			
Volume of HCl added in cm^3			

PRACTICAL MANUAL

B.Sc - IV Semester

2019-20

Name:.....

Seat No:.....



Department of Chemistry
SDM College
Ujire

CONTENTS

Sl. No.	Title	Page No.
1	Determination of rate constant of acid hydrolysis of methyl acetate	2- 5
2	Effect of catalytic strength of acids on hydrolysis of methyl acetate	6- 9
3	Effect of acid strength on hydrolysis of methyl acetate	10 - 13
4	Determination of molecular weight by Walker Lumsden method	14 - 17
5	Percentage of sodium chloride by CST method	18 - 21
6	Percentage of toluene by refractometer	22 - 25
7	Determination of density and surface tension of the given liquid	26 - 29
8	Determination of density and viscosity of the given organic liquid	30 - 33
9	Composition of glycerol water mixture	34 - 37
10	Catalytic decomposition of hydrogen peroxide	38 - 39
11	Verification of Distribution law	40 - 43
12	Determination of partition coefficient	44 - 47
13	Preparation of colloidal arsenious sulphide solution	49

Observations and Tabulation:

Time in min	Burette Readings		Volume of NaOH in cm ³	$k = \frac{2.303}{t} \log \frac{V_{\infty} - V_0}{V_{\infty} - V_t} \text{ min}^{-1}$
	Final	Initial		
0				
10				
20				
30				
40				
∞				

Mean value of $k = \dots \dots \dots \text{ min}^{-1}$

Hence, the rate constant of the reaction is = $\dots \dots \dots \text{ min}^{-1}$


Expt No.:

Date:

DETERMINATION OF VELOCITY CONSTANT OF ACID HYDROLYSIS OF METHYL ACETATE

Aim: To determine the rate constant of acid hydrolysis of methyl acetate

Principle: Methyl acetate is an ester. It is readily hydrolysed by dilute hydrochloric acid into methyl alcohol and acetic acid.



Two molecules of reactants are involved in the reaction. Since the quantity of ester is large, its concentration practically remains constant. Hence, the above reaction is of pseudo first order. Its rate constant can be calculated using the expression.

$$k = \frac{2.303}{t} \log \frac{a}{a-x}$$

where a is the initial concentration and x is the decrease in concentration of the ester in a time t . The progress of the reaction is followed by titrating a known volume of the solution against standard sodium hydroxide solution at regular intervals of time. Let V_0 , V_t and V_∞ represent the titre values in the beginning, after a time t and at the end of the reaction respectively.

$V_t - V_0 \propto$ Amount of acetic acid formed in a time t (x)

$V_\infty - V_0 \propto$ Amount of acetic acid formed at the end (It is the same as the amount of ester taken in the beginning) (a)

$$(V_\infty - V_0) - V_t - V_0 \propto a - x$$

$$(V_\infty - V_t) \propto a - x$$

The above first order equation becomes,

$$k = \frac{2.303}{t} \log \frac{V_\infty - V_0}{V_\infty - V_t} \text{ min}^{-1}$$

Chemicals Required:

100 cm³ of 0.5 N HCl, 15 cm³ of methyl acetate and 0.1 N NaOH

Procedure: 100 cm³ of 0.5N HCl is taken in a dry 250 cm³ stoppered bottle. This bottle is kept in a water bath which acts as a thermostat. 15cm³ of methyl acetate taken in another stoppered bottle is also kept in the water bath

A burette is filled with 0.1N NaOH. When the solution in the bottles have attained constant temperature, 5cm³ of the ester is pipetted out into the bottle containing the hydrochloric acid and shaken well.

5cm³ of the reaction mixture is immediately pipetted out into a conical flask containing ice cold water (to arrest further reaction) and a drop of phenolphthalein is added. When the pipette is half empty a stop clock is started. The solution is titrated against 0.1N NaOH. This titre value gives the amount of HCl present in the reaction mixture (V_0).

The titrations are repeated with the same volume of the mixture after 10, 20, 30 and 40 minutes, in the same manner to get V_t .

After the first titration, a portion of (about 25 cm³) the reaction mixture is taken in a stoppered bottle and kept in a water bath at 60⁰ C for 1 hour to complete the reaction. It is cooled to the original temperature. 5cm³ of the solution is titrated against the same NaOH and cooled to the original temperature. 5 cm³ of the solution is titrated against the same NaOH solution to get V_∞

Result:

The value of k is found to be a constant. Therefore the reaction is of first order.

The velocity constant of acid hydrolysis of methyl acetate =min⁻¹

Observations and Tabulation

(a) Hydrolysis of methyl acetate using 0.5N HCl

Time in min	Burette Readings		Volume of NaOH in cm ³	$k_t = \frac{2.303}{t} \log \frac{V_\infty - V_0}{V_\infty - V_t} \text{ min}^{-1}$
	Final	Initial		
0				
10				
20				
30				
40				
∞				

Mean value of $k_t = \dots \dots \dots \text{ min}^{-1}$

Hence, the rate constant of the reaction using 0.5N HCl is

$= \dots \dots \dots \text{ min}^{-1}$



Sri Dharmasthala Manjunatheshwara College
(Autonomous), Ujire-574 240, Dakshina Kannada, Karnataka State

PRACTICAL MANUAL

B.Sc- VI Semester

2019-20

Name:.....

Seat No:.....




Department of Chemistry

SDM (Autonomous) College

Ujire

INDEX

Serial No.	Name of the experiment	Page No
1.	Acetanilide from aniline	2-3
2.	Benzoic acid from toluene	4-5
3.	Meta dinitrobenzene from nitrobenzene	6-7
4.	Para bromo acetanilide from acetanilide	8-9
5.	Preparation of 1,1-bis-2-naphthol	10-11
6.	P- bromo acetanilide from acetanilide (green Chemistry)	12-13
7.	Tetraammine copper(II)sulphate	14-15
8.	Potassium trioxalato ferrate(III)	16-17
9.	Hexamine cobalt (III) chloride	18-19
10.	Inversion of cane sugar	20-23
11.	Equivalent Conductance at Infinite Dilution	24-27
12.	Potentiometric Titration	28-33
13.	Colorimetry	34-37
14.	Conductometric Titration	38-41
15.	Dissociation Constant of Weak Acid	42-45
16.	Microscale preparation P - bromo acetanilide from acetanilide	46-47



Expt.No:

Date:

Acetanilide from Aniline

Aim: To prepare acetanilide from aniline.

Principle: Aniline is acetylated with acetic anhydride to get acetanilide.


Procedure: 2.5cm³ Aniline is taken in a clean and dry 100ml conical flask. Add about 8ml of pure acetic anhydride into it in a bulk with stirring.

The contents of the flask were swirled for about 15-20 min. The contents of the flask or the semi liquid of reaction mixture is poured into about 150ml of the ice water taken in a 400 ml beaker with constant stirring. The crude sample of acetanilide is filtered and washed with ice cold water until the washings are free from acetic acid. The crude product is dried and weighed.

A portion of the sample is re crystallized from hot water. The sample is dried between the folds of the filter paper and melting point is determined.

Result: Yield of the product =g

Melting point = °C



Toluene

Benzoic acid

Expt.No:

Date:

Benzoic acid from toluene

Aim: To prepare benzoic acid from toluene

Principle: Toluene is refluxed with alkaline potassium permanganate solution and the sodium salt obtained is acidified with dilute hydrochloric acid to get benzoic acid.


Procedure: 2.5cm^3 of toluene is taken in a round bottomed flask fitted with water cooled reflux condenser. 5cm^3 of alkaline potassium permanganate (1g of sodium carbonate + 2g of potassium permanganate is dissolved in 50cc of water) is added to the flask. 5cm^3 of above solution is added to the round bottomed flask at definite intervals of time, till the 50cm^3 of solution is exhausted.

During this time, till the contents of round bottomed flask is refluxed. When the addition of alkaline potassium permanganate is completed, it is found that the oily layer in the flask is disappeared. Now the flask is cooled and the reflux condenser is disconnected. The contents of the flask are filtered to remove the brown manganous dioxide. The filtrate is now acidified by adding 5cm^3 of conc. Hydrochloric acid. A little of sodium sulphite is added with stirring (if any brown colour is found in the solution), until the colour completely disappears. Now the white crystals of benzoic acid are separated from the solution by filtration. It is then washed with cold water and dried between the folds of filter paper. Yield of the crude product is determined.

A portion of crude product is recrystallised from hot water. The recrystallised sample is dried and the melting point is determined.

Result: Yield of the product =g

Melting point = $^{\circ}\text{C}$



ATTESTED


 Principal
 PRINCIPAL
 S.D.M COLLEGE
 UJIRE - 574240, D.K.